

For Research Use Only

RPA1 Monoclonal antibody, PBS Only

Catalog Number: 67973-1-PBS



Basic Information

Catalog Number:

67973-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28239

GenBank Accession Number:

BC018126

GeneID (NCBI):

6117

UNIPROT ID:

P27694

Full Name:

replication protein A1, 70kDa

Calculated MW:

616 aa, 68 kDa

Observed MW:

68 kDa

Purification Method:

Protein G purification

CloneNo.:

1G3F11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67973-1-Ig (RPA1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Beta Actin Monoclonal antibody (HRP-66009) as loading control. This data was developed using the same antibody clone with 67973-1-PBS in a different storage buffer formulation.