

For Research Use Only

KCNA1 Monoclonal antibody, PBS Only



Catalog Number: 68008-1-PBS

Basic Information

Catalog Number:

68008-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG12563

GenBank Accession Number:

BC101733

GeneID (NCBI):

3736

UNIPROT ID:

Q09470

Full Name:

potassium voltage-gated channel, shaker-related subfamily, member 1 (episodic ataxia with myokymia)

Calculated MW:

495 aa, 56 kDa

Observed MW:

75 kDa

Purification Method:

Protein A purification

CloneNo.:

1D6B9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Pig, Mouse, Rat, Rabbit

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68008-1-Ig (KCNA1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68008-1-PBS in a different storage buffer formulation.