

For Research Use Only

RALY Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68011-3-PBS



Basic Information

Catalog Number: 68011-3-PBS	GenBank Accession Number: BC103753	Purification Method: Protein A Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 22913	CloneNo.: 2H1H10
Source: Mouse	UNIPROT ID: Q9UKM9	
Isotype: IgG3	Full Name: RNA binding protein, autoantigenic (hnRNP-associated with lethal yellow homolog (mouse))	
Immunogen Catalog Number: AG12250	Calculated MW: 307 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

68011-3-PBS targets RALY as part of a matched antibody pair.

MP51385-1: 68011-2-PBS capture and 68011-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

