

For Research Use Only

Endophilin 1 Monoclonal antibody

Catalog Number: 68029-1-Ig



Basic Information

Catalog Number: 68029-1-Ig	GenBank Accession Number: BC032825	Purification Method: Protein A purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6456	CloneNo.: 1G4H1
Source: Mouse	Full Name: SH3-domain GRB2-like 2	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG2a	Calculated MW: 39.8 kDa	
Immunogen Catalog Number: AG24793	Observed MW: 38-40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: Human, mouse, rat, pig, rabbit	WB : pig brain tissue, rabbit brain, rat brain, mouse brain, rabbit cerebellum, rat cerebellum

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68029-1-Ig (SH3GL2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.