

For Research Use Only

Endophilin 1 Monoclonal antibody

Catalog Number: 68029-1-Ig



Basic Information

Catalog Number:

68029-1-Ig

Size:

150ul , Concentration: 1000 ug/ml by

Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG24793

GenBank Accession Number:

BC032825

GeneID (NCBI):

6456

UNIPROT ID:

Q99962

Full Name:

SH3-domain GRB2-like 2

Calculated MW:

39.8 kDa

Observed MW:

38-40 kDa

Purification Method:

Protein A purification

CloneNo.:

1G4H1

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat, pig, rabbit

Positive Controls:

WB : pig brain tissue, rabbit brain, rat brain, mouse brain, rabbit cerebellum, rat cerebellum

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68029-1-Ig (SH3GL2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.