

For Research Use Only

RAB12 Monoclonal antibody

Catalog Number: 68079-1-Ig **1 Publications**



Basic Information

Catalog Number: 68079-1-Ig	GenBank Accession Number: NM_001025300	Purification Method: Protein A purification
Size: 150ul , Concentration: 1000 ug/ml by Nanodrop;	GeneID (NCBI): 201475	CloneNo.: 2G11C2
Source: Mouse	UNIPROT ID: Q6IQ22	Recommended Dilutions: WB: 1:2000-1:10000
Isotype: IgG2a	Full Name: RAB12, member RAS oncogene family	
Immunogen Catalog Number: AG31627	Calculated MW: 27KD	
	Observed MW: 27-30 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HEK-293 cells, HepG2 cells, HUVEC cells, A549 cells, LNCaP cells, U-251 cells, HSC-T6 cells, RAW 264.7 cells
Cited Applications: WB	
Species Specificity: Human, mouse, rat	
Cited Species: human	

Notable Publications

Author	Pubmed ID	Journal	Application
Yuka Suimon	40204187	Am J Pathol	WB

Storage


Storage:
Store at -20°C.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68079-1-Ig (RAB12 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.