

For Research Use Only

RAB12 Monoclonal antibody

Catalog Number: 68079-1-Ig



Basic Information

Catalog Number:

68079-1-Ig

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG31627

GenBank Accession Number:

NM_001025300

GeneID (NCBI):

201475

UNIPROT ID:

Q6IQ22

Full Name:

RAB12, member RAS oncogene family

Calculated MW:

27KD

Observed MW:

27-30 kDa

Purification Method:

Protein A purification

CloneNo.:

2G11C2

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : HeLa cells, HEK-293 cells, HepG2 cells, HUVEC cells, A549 cells, LNCaP cells, U-251 cells, HSC-T6 cells, RAW 264.7 cells

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68079-1-Ig (RAB12 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.