

For Research Use Only

# PYCRL Monoclonal antibody, PBS Only



Catalog Number: 68087-1-PBS

## Basic Information

Catalog Number:

68087-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31938

GenBank Accession Number:

BC007993

GeneID (NCBI):

65263

UNIPROT ID:

Q53H96

Full Name:

pyrroline-5-carboxylate reductase-like

Calculated MW:

274 aa, 29 kDa

Observed MW:

26-29 kDa

Purification Method:

Protein G purification

CloneNo.:

1G7A3

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Rabbit

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)


W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68087-1-Ig (PYCRL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68087-1-PBS in a different storage buffer formulation.



rabbit kidney tissue were subjected to SDS PAGE followed by western blot with 68087-1-Ig (PYCRL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68087-1-PBS in a different storage buffer formulation.