

For Research Use Only

RPL12 Monoclonal antibody, PBS Only

Catalog Number: 68110-1-PBS



Basic Information

Catalog Number:

68110-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6580

GenBank Accession Number:

BC050644

GeneID (NCBI):

6136

UNIPROT ID:

P30050

Full Name:

ribosomal protein L12

Calculated MW:

18 kDa

Observed MW:

21 kDa

Purification Method:

Protein G purification

CloneNo.:

2C9B3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68110-1-Ig (RPL12 antibody) at dilution of 1:100000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68110-1-PBS in a different storage buffer formulation.