

For Research Use Only

P30 of ASFV Monoclonal antibody

Catalog Number: 68134-1-Ig



Basic Information

Catalog Number:

68134-1-Ig

GenBank Accession Number:

M96354

Purification Method:

Protein G purification

Size:

150ul , Concentration: 1000 ug/ml by
Nanodrop;

GeneID (NCBI):

22220322

CloneNo.:

3B9E3

Source:

Mouse

UNIPROT ID:

P34204

Recommended Dilutions:

WB 1:5000-1:50000

Isotype:

IgG1

Full Name:

Phosphoprotein p30

Calculated MW:

30 kDa

Immunogen Catalog Number:

AG31446

Observed MW:

35 kDa

Applications

Tested Applications:

WB, ELISA

Positive Controls:

WB : Ag31446,

Species Specificity:

Human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein of P30 of ASFV