

For Research Use Only

P30 of ASFV Monoclonal antibody

Catalog Number: 68134-1-Ig



Basic Information

Catalog Number: 68134-1-Ig	GenBank Accession Number: M96354	Purification Method: Protein G purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 22220322	CloneNo.: 3B9E3
Source: Mouse	Full Name: Phosphoprotein p30	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG1	Calculated MW: 30 kDa	
Immunogen Catalog Number: AG31446	Observed MW: 35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Ag31446,
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein of P30 of ASFV