

For Research Use Only

Asfv P30 of ASFV Monoclonal antibody, PBS Only

Catalog Number: 68134-1-PBS

Basic Information

Catalog Number:

68134-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG31446

GenBank Accession Number:

M96354

GeneID (NCBI):

22220322

UNIPROT ID:

P34204

Full Name:

Phosphoprotein p30

Calculated MW:

30 kDa

Observed MW:

35 kDa

Purification Method:

Protein G purification

CloneNo.:

3B9E3

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein of P30 of ASFV This data was developed using the same antibody clone with 68134-1-PBS in a different storage buffer formulation.