

For Research Use Only

EPS8L2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68330-2-PBS

Basic Information

Catalog Number: 68330-2-PBS	GenBank Accession Number: BC002474	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 64787	CloneNo.: 3B10G8
Source: Mouse	UNIPROT ID: Q9H6S3	
Isotype: IgG1	Full Name: EPS8-Like 2	
Immunogen Catalog Number: AG16808	Calculated MW: 715 aa, 81 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

68330-2-PBS targets EPS8L2 as part of a matched antibody pair:

MP50248-1: 68330-2-PBS capture and 68330-3-PBS detection (validated in Cytometric bead array)

MP50248-2: 68330-2-PBS capture and 68330-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50248-1, EPS8L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68330-2-PBS. Detection antibody: 68330-3-PBS. Standard: Ag16808. Range: 0.781-200 ng/mL.



Cytometric bead array standard curve of MP50248-2, EPS8L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68330-2-PBS. Detection antibody: 68330-4-PBS. Standard: Ag16808. Range: 0.391-200 ng/mL.