

For Research Use Only

# EPS8L2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68330-3-PBS

## Basic Information

Catalog Number: 68330-3-PBS	GenBank Accession Number: BC002474	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 64787	CloneNo.: 2G5E6
Source: Mouse	UNIPROT ID: Q9H6S3	
Isotype: IgG1	Full Name: EPS8-Like 2	
Immunogen Catalog Number: AG16808	Calculated MW: 715 aa, 81 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Product Information

68330-3-PBS targets EPS8L2 as part of a matched antibody pair:

MP50248-1: 68330-2-PBS capture and 68330-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50248-1, EPS8L2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68330-2-PBS. Detection antibody: 68330-3-PBS. Standard:Ag16808. Range: 0.781-200 ng/mL.