

For Research Use Only

MITD1 Monoclonal antibody, PBS Only



Catalog Number: 68367-1-PBS

Basic Information

Catalog Number:

68367-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11237

GenBank Accession Number:

BC018453

GeneID (NCBI):

129531

UNIPROT ID:

Q8WV92

Full Name:

MIT, microtubule interacting and transport, domain containing 1

Calculated MW:

249 aa, 29 kDa

Observed MW:

30 kDa

Purification Method:

Protein G purification

CloneNo.:

1F11E2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, Mouse, Rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68367-1-Ig (MITD1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Beta Actin Monoclonal antibody (HRP-66009) as loading control. This data was developed using the same antibody clone with 68367-1-PBS in a different storage buffer formulation.