

For Research Use Only

ZNF189 Monoclonal antibody, PBS Only



Catalog Number: 68424-1-PBS

Basic Information

Catalog Number: 68424-1-PBS	GenBank Accession Number: BC126120	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7743	CloneNo.: 2A4B4
Source: Mouse	UNIPROT ID: O75820	
Isotype: IgG2b	Full Name: zinc finger protein 189	
Immunogen Catalog Number: AG33164	Observed MW: 70-73 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68424-1-Ig (ZNF189 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68424-1-PBS in a different storage buffer formulation.



HeLa cells were subjected to SDS PAGE followed by western blot with 68424-1-Ig (ZNF189 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68424-1-PBS in a different storage buffer formulation.