

For Research Use Only

ATP6V1A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68440-4-PBS



Basic Information

Catalog Number:

68440-4-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG10820

GenBank Accession Number:

BC013138

GeneID (NCBI):

523

UNIPROT ID:

P38606

Full Name:

ATPase, H⁺ transporting, lysosomal 70kDa, V1 subunit A

Calculated MW:

617 aa, 68 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2F11E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Product Information

68440-4-PBS targets ATP6V1A as part of a matched antibody pair:

MP51785-2: 68440-2-PBS capture and 68440-4-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51785-2, ATP6V1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68440-2-PBS. Detection antibody: 68440-4-PBS.



Cytometric bead array standard curve of MP51785-2, ATP6V1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68440-2-PBS. Detection antibody: 68440-4-PBS. Standard: Ag10820. Range: 0.781-100 ng/mL