

For Research Use Only

# EBP2 Monoclonal antibody, PBS Only

Catalog Number: 68477-1-PBS



## Basic Information

Catalog Number:

68477-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8833

GenBank Accession Number:

BC009175

GeneID (NCBI):

10969

UNIPROT ID:

Q99848

Full Name:

EBNA1 binding protein 2

Calculated MW:

306 aa, 35 kDa

Observed MW:

38-40 kDa

Purification Method:

Protein A purification

CloneNo.:

1F6G1

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68477-1-Ig (EBP2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68477-1-PBS in a different storage buffer formulation.