

For Research Use Only

# DIABLO Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68480-4-PBS

## Basic Information

<b>Catalog Number:</b> 68480-4-PBS	<b>GenBank Accession Number:</b> BC011909	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 56616	<b>CloneNo.:</b> 2H1A1
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9NR28	
<b>Isotype:</b> IgG2b	<b>Full Name:</b> diablo homolog (Drosophila)	
<b>Immunogen Catalog Number:</b> AG32774	<b>Calculated MW:</b> 27 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

68480-4-PBS targets DIABLO as part of a matched antibody pair:

MP50334-2: 68480-4-PBS capture and 68480-1-PBS detection (validated in Cytometric bead array)

MP50334-3: 68480-4-PBS capture and 68480-5-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50334-2, DIABLO Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68480-4-PBS. Detection antibody: 68480-1-PBS. Standard:Ag32774. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50334-3, DIABLO Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68480-4-PBS. Detection antibody: 68480-5-PBS. Standard:Ag32774. Range: 0.781-50 ng/mL.