

For Research Use Only

FLAD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68491-7-PBS



Basic Information

Catalog Number: 68491-7-PBS	GenBank Accession Number: BC011378	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 80308	CloneNo.: 1F1H1
Source: Mouse	UNIPROT ID: Q8NFF5	
Isotype: IgG1	Full Name: FAD1 flavin adenine dinucleotide synthetase homolog (S. cerevisiae)	
Immunogen Catalog Number: AG6845	Calculated MW: 446 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68491-7-PBS targets FLAD1 as part of a matched antibody pair.

MP50710-4: 68491-7-PBS capture and 68491-6-PBS detection (validated in Cytometric bead array)

MP50710-5: 68491-7-PBS capture and 68491-1-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50710-4, FLAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68491-7-PBS. Detection antibody: 68491-6-PBS. Standard: Ag6845. Range: 6.25-100 ng/mL



Cytometric bead array standard curve of MP50710-5, FLAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68491-7-PBS. Detection antibody: 68491-1-PBS. Standard: Ag6845. Range: 1.563-100 ng/mL