

For Research Use Only

C7orf50 Monoclonal antibody, PBS Only



Catalog Number: 68494-1-PBS

Basic Information

Catalog Number: 68494-1-PBS	GenBank Accession Number: BC006224	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 84310	CloneNo.: 2A3G5
Source: Mouse	UNIPROT ID: Q9BRJ6	
Isotype: IgG2a	Full Name: chromosome 7 open reading frame 50	
Immunogen Catalog Number: AG15142	Calculated MW: 194 aa, 22 kDa	
	Observed MW: 22 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

C7orf50 (Chromosome 7, Open Reading Frame 50) is ubiquitously expressed in the kidneys, brain, fat, prostate, spleen, among 22 other tissues and demonstrates low tissue specificity.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68494-1-Ig (C7orf50 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68494-1-PBS in a different storage buffer formulation.