

For Research Use Only

RAD17 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68539-2-PBS



Basic Information

| | | |
|---|---|--|
| Catalog Number: 68539-2-PBS | GenBank Accession Number: BC032304 | Purification Method: Protein G Magarose purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 5884 | CloneNo.: 2G8G5 |
| Source: Mouse | UNIPROT ID: O75943 | |
| Isotype: IgG1 | Full Name: RAD17 homolog (S. pombe) | |
| Immunogen Catalog Number: AG4348 | Calculated MW: 681 aa, 71 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68539-2-PBS targets RAD17 as part of a matched antibody pair:

MP50647-1: 68539-2-PBS capture and 68539-3-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50647-1, RAD17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68539-2-PBS. Detection antibody: 68539-3-PBS. Standard:Ag4348. Range: 0.391-12.5 ng/mL