

For Research Use Only

MRPL24 Monoclonal antibody, PBS Only



Catalog Number: 68556-1-PBS

Basic Information

Catalog Number:

68556-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9576

GenBank Accession Number:

BC016700

GeneID (NCBI):

79590

UNIPROT ID:

Q96A35

Full Name:

mitochondrial ribosomal protein L24

Calculated MW:

216 aa, 25 kDa

Observed MW:

25-28 kDa

Purification Method:

Protein G purification

CloneNo.:

3D3C8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68556-1-Ig (MRPL24 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68556-1-PBS in a different storage buffer formulation.