

For Research Use Only

TRPV2 Monoclonal antibody, PBS Only



Catalog Number: 68563-1-PBS

Basic Information

Catalog Number:

68563-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG8855

GenBank Accession Number:

BC018926

GeneID (NCBI):

51393

UNIPROT ID:

Q9Y5S1

Full Name:

transient receptor potential cation channel, subfamily V, member 2

Calculated MW:

764 aa, 86 kDa

Observed MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

2F11H3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68563-1-Ig (TRPV2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68563-1-PBS in a different storage buffer formulation.