

For Research Use Only

PIP5K3 Monoclonal antibody, PBS Only



Catalog Number: 68571-1-PBS

Basic Information

Catalog Number:

68571-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4366

GenBank Accession Number:

BC032389

GeneID (NCBI):

200576

UNIPROT ID:

Q9Y2I7

Full Name:

phosphatidylinositol-3-phosphate/phosphatidylinositol 5-kinase, type III

Calculated MW:

451 aa, 50 kDa, 237 kDa

Observed MW:

240-250 kDa

Purification Method:

Protein G purification

CloneNo.:

2F5G1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat, rabbit, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68571-1-Ig (PIP5K3 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68571-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 68571-1-Ig (PIP5K3 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68571-1-PBS in a different storage buffer formulation.