

For Research Use Only

TRADD Monoclonal antibody, PBS Only



Catalog Number: 68573-1-PBS

Basic Information

Catalog Number:

68573-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7768

GenBank Accession Number:

BC004491

GeneID (NCBI):

8717

UNIPROT ID:

Q15628

Full Name:

TNFRSF1A-associated via death domain

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Protein G purification

CloneNo.:

1C3H3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, pig, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68573-1-Ig (TRADD antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68573-1-PBS in a different storage buffer formulation.