

For Research Use Only

SSH3 Monoclonal antibody

Catalog Number: 68583-1-Ig



Basic Information

Catalog Number: 68583-1-Ig	GenBank Accession Number: BC004210	Purification Method: Protein A purification
Size: 150ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 54961	CloneNo.: 1A9C9
Source: Mouse	Full Name: slingshot homolog 3 (Drosophila)	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG2b	Calculated MW: 659 aa, 73 kDa	
Immunogen Catalog Number: AG13139	Observed MW: 90-95 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Positive Controls:

WB: A431 cells, A549 cells, T-47D cells, MCF-7 cells, HeLa cells, HT-29 cells, HepG2 cells, HuH-7 cells, THP-1 cells

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68583-1-Ig (SSH3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.