

For Research Use Only

SSH3 Monoclonal antibody, PBS Only

Catalog Number: 68583-1-PBS



Basic Information

Catalog Number:

68583-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG13139

GenBank Accession Number:

BC004210

GeneID (NCBI):

54961

UNIPROT ID:

Q8TE77

Full Name:

slingshot homolog 3 (Drosophila)

Calculated MW:

659 aa, 73 kDa

Observed MW:

90-95 kDa

Purification Method:

Protein A purification

CloneNo.:

1A9C9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68583-1-Ig (SSH3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68583-1-PBS in a different storage buffer formulation.