

For Research Use Only

TXNDC5 Monoclonal antibody, PBS Only

Catalog Number: 68609-1-PBS



Basic Information

Catalog Number:

68609-1-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13965

GenBank Accession Number:

BC001199

GeneID (NCBI):

81567

UNIPROT ID:

Q8NBS9

Full Name:

thioredoxin domain containing 5 (endoplasmic reticulum)

Calculated MW:

432 aa, 48 kDa

Observed MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

1A12C3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68609-1-Ig (TXNDC5 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68609-1-PBS in a different storage buffer formulation.