

For Research Use Only

MRPL10 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68623-2-PBS



Basic Information

Catalog Number: 68623-2-PBS	GenBank Accession Number: BC052601	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 124995	CloneNo.: 3C2F4
Source: Mouse	UNIPROT ID: Q7Z7H8	
Isotype: IgG1	Full Name: mitochondrial ribosomal protein L10	
Immunogen Catalog Number: AG10060	Calculated MW: 261 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68623-2-PBS targets MRPL10 as part of a matched antibody pair:

MP50459-1: 68623-2-PBS capture and 68623-3-PBS detection (validated in Cytometric bead array)

MP50459-2: 68623-2-PBS capture and 68623-4-PBS detection (validated in Cytometric bead array)

MP50459-4: 68623-5-PBS capture and 68623-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50459-1, MRPL10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68623-2-PBS. Detection antibody: 68623-3-PBS. Standard:Ag10060. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50459-2, MRPL10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68623-2-PBS. Detection antibody: 68623-4-PBS. Standard:Ag10060. Range: 0.098-50 ng/mL



Cytometric bead array standard curve of MP50459-4, MRPL10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68623-5-PBS. Detection antibody: 68623-2-PBS. Standard:Ag10060. Range: 0.098-100 ng/mL