

For Research Use Only

DPF3 Monoclonal antibody, PBS Only

Catalog Number: 68646-1-PBS



Basic Information

Catalog Number:

68646-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG33877

GenBank Accession Number:

BC060801

GeneID (NCBI):

8110

UNIPROT ID:

Q92784

Full Name:

D4, zinc and double PHD fingers, family 3

Calculated MW:

357 aa, 40 kDa

Observed MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

1F3H4

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, rabbit, rat, mouse, pig

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68646-1-Ig (DPF3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68646-1-PBS in a different storage buffer formulation.