

For Research Use Only

# CD200 Monoclonal antibody, PBS Only

Catalog Number: 68700-3-PBS



## Basic Information

<b>Catalog Number:</b> 68700-3-PBS	<b>GenBank Accession Number:</b> BC022522	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 4345	<b>CloneNo.:</b> 2G8D1
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P41217	
<b>Isotype:</b> IgG1	<b>Full Name:</b> CD200 molecule	
<b>Immunogen Catalog Number:</b> EG0006	<b>Calculated MW:</b> 269 aa, 30 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA

**Species Specificity:**  
human

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data