

For Research Use Only

# LDB1 Monoclonal antibody, PBS Only

Catalog Number: 68721-1-PBS



## Basic Information

Catalog Number:

68721-1-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34563

GenBank Accession Number:

BC000482

GeneID (NCBI):

8861

UNIPROT ID:

Q86U70

Full Name:

LIM domain binding 1

Calculated MW:

46 kDa

Purification Method:

Protein G purification

CloneNo.:

2F11C3

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68721-1-Ig (LDB1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68721-1-PBS in a different storage buffer formulation.