

For Research Use Only

Medin Monoclonal antibody

Catalog Number: 68736-1-Ig



Basic Information

Catalog Number:

68736-1-Ig

Size:

150ul , Concentration: 1000 ug/ml by

Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30145

GenBank Accession Number:

NM_005928

GeneID (NCBI):

4240

UNIPROT ID:

Q08431

Full Name:

milk fat globule-EGF factor 8 protein

Calculated MW:

43 kDa

Observed MW:

51 kDa

Purification Method:

Protein G purification

CloneNo.:

1H3D7

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human placenta tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68736-1-Ig (MFGE8 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.