

For Research Use Only

MFG E8 Monoclonal antibody, PBS Only



Catalog Number: 68736-1-PBS

Basic Information

Catalog Number: 68736-1-PBS	GenBank Accession Number: NM_005928	Purification Method: Protein G purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 4240	CloneNo.: 1H3D7
Source: Mouse	UNIPROT ID: Q08431	
Isotype: IgG1	Full Name: milk fat globule-EGF factor 8 protein	
Immunogen Catalog Number: AG30145	Calculated MW: 43 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68736-1-Ig (MFGE8 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68736-1-PBS in a different storage buffer formulation.