

For Research Use Only

IL-1 alpha Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68758-1-PBS



Basic Information

Catalog Number: 68758-1-PBS	GenBank Accession Number: BC013142	Purification Method: Protein G purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 3552	CloneNo.: 1B9F5
Source: Mouse	UNIPROT ID: P01583	
Isotype: IgG1	Full Name: interleukin 1, alpha	
Immunogen Catalog Number: HZ-1320	Calculated MW: 271 aa, 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68758-1-PBS targets IL-1 alpha as part of a matched antibody pair:

MP50113-1: 68757-1-PBS capture and 68758-1-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50113-1, IL1A/IL-1 alpha Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68757-1-PBS. Detection antibody: 68758-1-PBS. Standard:HZ-1320. Range: 0.098-3.125 ng/mL