

For Research Use Only

PPP2R1A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68778-1-PBS

Basic Information

Catalog Number: 68778-1-PBS	GenBank Accession Number: BC001537	Purification Method: Protein G Magarose purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 5518	CloneNo.: 1D1C1
Source: Mouse	UNIPROT ID: P30153	
Isotype: IgG1	Full Name: protein phosphatase 2 (formerly 2A), regulatory subunit A, alpha isoform	
Immunogen Catalog Number: AG8661	Calculated MW: 589 aa, 65 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

68778-1-PBS targets PPP2R1A as part of a matched antibody pair:

MP50144-1: 68778-1-PBS capture and 68778-2-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50144-1, PPP2R1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68778-1-PBS. Detection antibody: 68778-2-PBS. Standard:Ag8661. Range: 6.25-100 ng/mL