

For Research Use Only

Histone H3.1 Monoclonal antibody, PBS Only



Catalog Number: 68825-1-PBS

Basic Information

Catalog Number:

68825-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC066245

GeneID (NCBI):

8350

UNIPROT ID:

P68431

Full Name:

histone cluster 1, H3a

Purification Method:

Protein A purification

CloneNo.:

2E6D11

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68825-1-Ig (Histone H3.1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68825-1-PBS in a different storage buffer formulation.