

For Research Use Only

# IGFBP4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68834-3-PBS

## Basic Information

<b>Catalog Number:</b> 68834-3-PBS	<b>GenBank Accession Number:</b> BC016041	<b>Purification Method:</b> Protein G Magarose purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 3487	<b>CloneNo.:</b> 1G3C5
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P22692	
<b>Isotype:</b> IgG1	<b>Full Name:</b> insulin-like growth factor binding protein 4	
<b>Immunogen Catalog Number:</b> AG34800	<b>Calculated MW:</b> 28 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Product Information

68834-3-PBS targets IGFBP4 as part of a matched antibody pair:

MP50212-2: 68834-3-PBS capture and 68834-4-PBS detection (validated in Sandwich ELISA)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50212-2, IGFBP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68834-3-PBS. Detection antibody: HRP-conjugated 68834-4-PBS. Standard: Ag34800. Range: 0.049-50 ng/mL.