

For Research Use Only

KIF2A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68917-4-PBS

Basic Information

Catalog Number:

68917-4-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34099

GenBank Accession Number:

BC031828

GeneID (NCBI):

3796

UNIPROT ID:

O00139

Full Name:

kinesin heavy chain member 2A

Calculated MW:

679 aa, 77 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C3C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50342-3, KIF2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68917-1-PBS. Detection antibody: 68917-4-PBS. Standard:Ag34099. Range: 3.125-200 ng/mL.