

For Research Use Only

DNMT3B Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68919-2-PBS

Basic Information

Catalog Number:

68919-2-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33907

GenBank Accession Number:

NM_001207055

GeneID (NCBI):

1789

UNIPROT ID:

Q9UBC3

Full Name:

DNA (cytosine-5-)-methyltransferase 3 beta

Calculated MW:

95 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A3C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50346-1, DNMT3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68919-1-PBS. Detection antibody: 68919-2-PBS. Standard:Ag33907. Range: 0.391-200 ng/mL.