

For Research Use Only

CKAP4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68938-2-PBS

Basic Information

| | | |
|---|--|--|
| Catalog Number: 68938-2-PBS | GenBank Accession Number: BC082972 | Purification Method: Protein G Magarose purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 10970 | CloneNo.: 2B3H3 |
| Source: Mouse | UNIPROT ID: Q07065 | |
| Isotype: IgG1 | Full Name: cytoskeleton-associated protein 4 | |
| Immunogen Catalog Number: AG10022 | Calculated MW: 602 aa, 66 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

68938-2-PBS targets CKAP4 as part of a matched antibody pair:

MP50368-1: 68938-1-PBS capture and 68938-2-PBS detection (validated in Cytometric bead array)

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50368-1, CKAP4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68938-1-PBS. Detection antibody: 68938-2-PBS. Standard:Ag10022. Range: 1.563-100 ng/mL.