

For Research Use Only

IFN alpha 2B Monoclonal antibody

Catalog Number: 69016-1-Ig



Basic Information

Catalog Number: 69016-1-Ig	GenBank Accession Number: GeneID (NCBI): 3440	Purification Method: Protein G purification
Size: 100ug	Full Name: interferon, alpha 2	CloneNo.: 1E5C4
Source: Mouse		
Isotype: IgG1		
Immunogen Catalog Number: HZ-1072		

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Lyophilized antibodies are stable for 1 year from the date of receipt if stored between (-20°C) and (-80°C). Upon reconstitution we recommend that the solution can be stored at(4°C) for short term or at(-20°C) to (-80°C) for long term. Repeated freeze thaw cycles should be avoided with reconstituted products.

Storage Buffer:
Sterile PBS.

Aliquoting is unnecessary for -20° C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data