

For Research Use Only

# NeutraKine® IL-17A Monoclonal antibody



Catalog Number: 69021-2-Ig

## Basic Information

Catalog Number:

69021-2-Ig

Size:

100ug

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1113

GenBank Accession Number:

GeneID (NCBI):

3605

ENSEMBL Gene ID:

ENSG00000112115

Full Name:

interleukin 17A

Purification Method:

Protein G purification

CloneNo.:

1E3H1

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Storage

Storage:

Lyophilized antibodies are stable for 1 year from the date of receipt if stored between (-20°C) and (-80°C). Upon reconstitution we recommend that the solution can be stored at(4°C) for short term or at(-20°C) to (-80°C) for long term. Repeated freeze thaw cycles should be avoided with reconstituted products.

Storage Buffer:

Sterile PBS.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

