

For Research Use Only

rat anti SNAP/CLIP-tag® Monoclonal antibody (6F9)



www.ptglab.com

Catalog Number: 6f9

2 Publications

Basic Information

Catalog Number:
6f9

Clone No.:
6F9

Applications:
WB

Host/IsoType:
Rat/IgG2a

Conjugate:
Unconjugated

Type:
Primary Antibody

Class:
Monoclonal

RRID:
AB_2827567

Purification Method:
Protein A or G purified

Description

Rat monoclonal antibody [6F9] to SNAP-tag and CLIP-tag

Specificity/Target

SNAP-tag- and CLIP-tagged fusion proteins

Recommended Dilution

Western Blot: 1:1,000

Storage

Storage:

Shipped on blue ice. Upon receipt store at -20°C. Stable for one year.

Storage Buffer:


PBS, 50% glycerol, 0.09% sodium azide

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Western blot analysis of cell extract from HEK293T cells transiently expressing SNAP-Tubulin and from untransfected cells.



Sensitivity tested with purified recombinant SNAP-tag protein