

For Research Use Only

rat anti GST Monoclonal antibody (6G9)



www.ptglab.com

Catalog Number: 6g9

4 Publications

Basic Information

Catalog Number:
6g9

Clone No.:
6G9

Applications:
Western Blot, ELISA

Host/IsoType:
Rat/IgG2a

Conjugate:
Unconjugated

Type:
Primary Antibody

Class:
Monoclonal

RRID:
AB_2631397

Purification Method:
Protein G purified

Description

Rat monoclonal [6G9] to Glutathione S-Transferase (GST) of Schistosoma japonicum

Specificity/Target

The antibody recognizes GST of Schistosoma japonicum, not tested on other species.

Recommended Dilution

Western Blot: 1:1,000
ELISA: 1:1,000

Storage


Storage:
Shipped on blue ice. Upon receipt store at 4°C; stable for one year
Storage Buffer:
PBS, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Western blot analysis of cell extract from BL21 cells expressing GST



Western blot, dilution series of purified recombinant GST-protein, anti-GST antibody