

For Research Use Only

rat anti HA Monoclonal antibody (7C9)



www.ptglab.com

Catalog Number: 7c9

28 Publications

Basic Information

Catalog Number:
7c9

Clone No.:
7C9

Applications:
WB, IF, IP

Host/IsoType:
Rat/IgG2a

Conjugate:
Unconjugated

Type:
Primary Antibody

Class:
Monoclonal

RRID:
AB_2631399

Purification Method:
Protein G purified

Description

Rat monoclonal [7C9] to HA-tag

Specificity/Target

HA-tag (YPYDVPDYA) at the N-terminus, C-terminus, or internal site of fusion protein.

Recommended Dilution

Western Blot: 1:1,000
Immunofluorescence: 1:1,000
Immunoprecipitation: 5 µg/ IP

Storage

Storage:
Shipped on ice. Store at -20°C/-4°F. Stable for 1 year after shipment.


Storage Buffer:
PBS, preservative: 0.09 % Sodium azide

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Western blot analysis of cell extract from HEK293T cells expressing GFP-HA (29.9 kDa)



Immunofluorescence of HeLa cells transiently expressing HA-TOMM70