

For Research Use Only

# Phospho-p70 S6 kinase (Thr389) Recombinant monoclonal antibody

Catalog Number: 80021-1-RR

1 Publications



## Basic Information

Catalog Number:

80021-1-RR

Size:

100ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_003161

GeneID (NCBI):

6198

UNIPROT ID:

P23443

Full Name:

ribosomal protein S6 kinase, 70kDa,  
polypeptide 1

Purification Method:

Protein A purification

CloneNo.:

4C8

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Xue Jiao	40367705	Vet Microbiol	WB

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

