

For Research Use Only

STAT2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 80061-4-PBS



Basic Information

Catalog Number: 80061-4-PBS	GenBank Accession Number: BC051284	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6773	CloneNo.: 242685E6
Source: Rabbit	UNIPROT ID: P52630	
Isotype: IgG	Full Name: signal transducer and activator of transcription 2, 113kDa	
Immunogen Catalog Number: AG10168	Calculated MW: 851 aa, 98 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

80061-4-PBS targets STAT2 as part of a matched antibody pair.

MP01853-1: 80061-3-PBS capture and 80061-4-PBS detection (validated in Cytometric bead array)

MP01853-2: 80061-2-PBS capture and 80061-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01853-1, STAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80061-3-PBS. Detection antibody: 80061-4-PBS. Standard: Ag10168. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01853-2, STAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80061-2-PBS. Detection antibody: 80061-4-PBS. Standard: Ag10168. Range: 1.563-200 ng/mL.