

For Research Use Only

TTYH1 Recombinant antibody, PBS Only



Catalog Number: 80263-1-PBS

Basic Information

Catalog Number: 80263-1-PBS	GenBank Accession Number: BC011347	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 57348	CloneNo.: 8G21
Source: Rabbit	UNIPROT ID: Q9H313	
Isotype: IgG	Full Name: tweety homolog 1 (Drosophila)	
Immunogen Catalog Number: AG25290	Calculated MW: 49-50 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, rabbit, pig

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80263-1-RR (TTYH1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80263-1-PBS in a different storage buffer formulation.